

Building reusable bioimaging workflows with Nextflow

Sebastian Gonzalez Tirado

Shout-outs

- James Fellow Yates (Leibniz-HKI / MPI-EVA)
- Christian Tischer (EMBL)

Outcomes of this talk

- What are workflow managers and why should I care?
- Examples of complex bioimaging workflows orchestrated with Nextflow.
- Practical recommendations to create your own workflows.

Have you ever faced one of these problems?

Discovering an exciting new tool

S: “Oh cool, let’s go to the repo and check it out!”



Cell segmentation in imaging-based spatial transcriptomics

Viktor Petukhov^{1,2}, Rosalind J. Xu^{1,2,3}, Ruslan A. Soldatov¹, Paolo Cadini^{1,2}, Konstantin Khodosevich^{1,2}, Jeffrey B. Moffitt^{1,2,3} and Peter V. Kharchenko^{1,2,3}

...could not even install it.

nature
biotechnology

ARTICLES

<https://doi.org/10.1038/s41587-020-0044-4>

Open Access

Repetitive use of tools



track
mate

Make sure your work is reusable

h-index from 0 to hero

Adapting slowly into standards

nf-core



Only one rule for bioimage analysts club: be lazy

What are workflow managers?

Like orchestra directors managing the flow of your pipeline/symphony.

What do they do?

Workflow managers provide a framework for the creation, execution, and monitoring of a pipeline. <...>

They simplify pipeline development, optimize resource usage, handle software installation and versions, and run on different compute platforms, enabling workflow portability and sharing.

Wratten et al. (2021) Nature Methods

Other benefits


You also benefit from

REUSABLE

- Increased **portability**
- More **efficiency**:
 - Less headaches from installation
 - Better resource management and monitoring
 - Extremely parallelised


Getting more technical

Nextflow is a reactive workflow framework and a programming DSL



Getting more technical

Nextflow is a reactive workflow framework and a programming DSL



In line with our first rule: do not reinvent the wheel


nf-core



A **community** effort to collect a curated set of analysis pipelines built using  Nextflow.

Introduction

nf-core/molkart is a pipeline for processing Molecular Cartography data from Resolve Bioscience (combinatorial FISH). It takes as input a table of FISH spot positions (x,y,z,gene), a corresponding DAPI image (TIFF format) and optionally an additional staining image in the TIFF format. nf-core/molkart performs end-to-end processing of the data including image processing, QC filtering of spots, cell segmentation, spot-to-cell assignment and reports quality metrics such as the spot assignment rate, average spots per cell and segmentation mask size ranges.



>_ run with

nf-core pipelines launch nf-cor

nf-core	Nextflow	Sequera Platform
---------	----------	------------------

subscribers	stars
163	13
open issues	open PRs
11	2
last release	last update
8 months ago	8 months ago

included modules

cellpose deepcell_mesmer ilastik_multicut
ilastik_pixelclassification
mindagap_duplicatefinder and 3 more modules

included subworkflows

utils_nextflow_pipeline utils_nfcore_pipeline
utils_nfschema_plugin

contributors



Adapting to nf-core specs takes time

You can also just write your own pipeline:

<https://github.com/BiolImageTools/iss-nf>

iss-nf: A Nextflow-based end-to-end *in situ* sequencing decoding workflow

 Nima Vakili, Sebastián González-Tirado,  Nils Kurzawa,  Dmytro Dvornikov, Zeinab Mokhtari, Frank Wippich,  Giovanna Bergamini,  Rainer Pepperkok,  Christian Tischer,  Luis A. Vale-Silva

doi: <https://doi.org/10.1101/2025.10.16.682795>

This article is a preprint and has not been certified by peer review [what does this mean?].

Main components of a nextflow script

```
images = Channel.fromPath("${params.inputDir}/*.tif")
    .map { file -> tuple(file.baseName, file) }

workflow = {
    labels = SEGMENT ( images )
    MEASURE ( images.join( labels ) )
}
```

```
workflow = {
    labels = SEGMENT ( images )
    MEASURE ( images.join( labels ) )
}
```

```
process SEGMENT {
    input:
    path(image)
    output:
    path("*.tif")

    script:
    """
    cellpose --input ${image} --pretrained_model cyto --save_tif
    """
}

process MEASURE {
    publishDir '${params.outputDir}', mode: 'copy'
    input:
    tuple path(image), path(labels)
    output:
    path("*.csv")

    script:
    """
    fiji --run measure.groovy "${image}, ${labels}, ${imageID}.csv"
    """
}
```

How do we run it?

\$ nextflow run your_nextflow_script.nf --inputDir /path/2/images

```
> nextflow run practice_nextflow.nf --inputDir /Users/sebgoti/Documents/PhD/Riken_Globias/training_repo/train_cyto2
Nextflow 25.04.8 is available - Please consider updating your version to it
```

```
N E X T F L O W ~ version 24.10.4
```

```
Launching `practice_nextflow.nf` [curious_lavoisier] DSL2 - revision: 4063a8a561
```

```
executor > local (10)
[1b/91e239] SEGMENT (10) [100%] 10 of 10 ✓
Completed at: 18-Oct-2025 08:48:55
Duration      : 16m 16s
CPU hours     : 1.6
Succeeded     : 10
```

Some of the most important aspects to consider

Configuration files (my_file.config):

```
process {  
  executor = 'slurm'  
  queue = 'htc-el8'  
  cpus = 8  
  memory = 16.GB  
  time = 20.min  
  container = "docker://segonzal/fish_analysis:0.0.3"  
  
  withLabel: 'long' {  
    cpus = 8  
    memory = 32.GB  
    container = "docker://segonzal/fish_analysis:0.0.3"  
  }  
  
  withLabel: 'registration' {  
    cpus = 4  
    memory = 32.GB  
  }  
  
  withLabel: 'min' {  
    cpus = 1  
    memory = 400.MB  
  }  
}
```

That is it for the theory!

Before moving on to the exercises, questions???

Familiarize yourself with it

1. Create your first nextflow script 'practice_nextflow.nf', define the channel with the input images and print this channel (tip: check the entry for the `.view()` operator on the [**nextflow documentation**](#)).

What you should see:

```
> nextflow run practice_nextflow.nf --inputDir /Users/sebgoti/Documents/PhD/Riken_Globias/training_repo/train_cyto2
Nextflow 25.04.8 is available - Please consider updating your version to it

NEXTFLOW ~ version 24.10.4

Launching `practice_nextflow.nf` [reverent_wiles] DSL2 - revision: 8e3c76088e

[649_img, /Users/sebgoti/Documents/PhD/Riken_Globias/training_repo/train_cyto2/649_img.png]
[731_img, /Users/sebgoti/Documents/PhD/Riken_Globias/training_repo/train_cyto2/731_img.png]
[659_img, /Users/sebgoti/Documents/PhD/Riken_Globias/training_repo/train_cyto2/659_img.png]
[721_img, /Users/sebgoti/Documents/PhD/Riken_Globias/training_repo/train_cyto2/721_img.png]
[596_img, /Users/sebgoti/Documents/PhD/Riken_Globias/training_repo/train_cyto2/596_img.png]
```


Familiarize yourself with it

2. Make your first process: use cellpose to segment the images from the channel you created in 1.

Familiarize yourself with it

3. Quantify the average intensity of the cells in the masks. You can create your own script or macro!